



TSC LIFE

Innovation with purpose

Cystoflex comes in two versions:

Cystoflex Standard Deflection	11010011
Cystoflex Reverse Deflection	11010012
Screeni (with/without accessories)	30030000 / 30030001

Note: Always consult the latest User Manual of Cystoflex devices.

1. Cystoflex User Manual ref. 11010144
2. Screeni User Manual ref. 30030511
3. 2024 Study: The Impact of Single-Use Cystoscopes on Clinical Time Workflow in an Outpatient Setting, (10.1016/j.urology.2024.03.033)

The use of our medical equipment is solely intended for the use of healthcare professionals. Healthcare professionals have the obligation to exercise their own professional clinical judgment when determining the suitability of a specific product for treating an individual patient. TSC Life does not provide medical advice and strongly advises healthcare professionals to undergo training on the usage of any specific product before using it. Healthcare professionals should consistently refer to the package insert, product label, and/or instructions for use, including guidelines for cleaning and sterilization (if applicable), before utilizing any product supplied by TSC Life. The information presented aims to display specific products and the extensive range of offerings from TSC Life. Availability of products may vary across markets due to regulatory and/or medical practices specific to each region. For inquiries regarding the availability of specific TSC Life products in your area, please reach out to your TSC Life representative.



TSC Life

Manufactured by : Axess Vision Technology, S.A.S. | Zone de La Liodière | 6 rue de la Flotière | 37300 Joué-lès-Tours | France
T +33 2 47343290 | letsconnect@tsc-life.com | www.tsc-life.com

TSC LIFE

CYSTOFLEX®



**Efficiency Without Compromise.
Smarter Cystoscopy, Seamless Workflow.**

tsc-life.com

CYSTOFLEX®

The Cystoflex single-use flexible cystoscope is designed to maintain the familiar feel of a reusable cystoscope while providing the added flexibility and convenience of a single-use solution: optimized for performance with a low environmental footprint.

Enhancing efficiencies



Ready-to-use solution
No need to worry about availability, Cystoflex is sterile and there for you when you need it.



Enhance patient safety
Cystoflex eliminates the risk of cross-contamination.



No reprocessing
Using Cystoflex frees up resources and streamlines workflow.



Cost-effective
No capital investment or repair costs.



Low Environmental Footprint : Made in France and weighing just 108g to minimize waste, Cystoflex is your best European option for reducing environmental impact and transportation emissions.

Also available from TSC Life: Fluido® Irrigation



During your endoscopic procedures, Fluido Irrigation is the worry-free and intuitive warming solution that empowers you to **optimize patient outcome and minimize risk of hypothermia.**

A controlled and stable temperature with no compromises on patient's safety.


Did you know?

+75%

Boost your efficiency with single-use cystoscopes and see up to 75% more patients per day compared to reusable scopes³.


Integrated T-Connector

Included for seamless instrument handling and a more efficient workflow. No hidden add-ons!




Intuitive Navigation

210° precise articulation that enables retroflexion to access challenging anatomies with ease.



Atraumatic Insertion Tip

Only 4mm, tapered and streamlined for smooth, painless introduction.



Exceptional Image Quality

Optimized brightness and full-screen mode ensure clear, detailed visuals with outstanding color rendering.